



recomLine Parvovirus B19 IgG [Avidity] recomLine Parvovirus B19 IgM

Strip-Immunoassay with antigens produced by recombinant techniques for the detection of IgG and IgM antibodies against human Parvovirus B19


Human Parvovirus B19 occurs endemically world-wide. At intervals of several years, it causes epidemics in small areas from late fall to early summer, affecting kindergartens and schools in particular.

Its surface is formed by a protein envelope which consists of two polypeptides (VP-1 and VP-2). VP-1 differs from VP-2 only by an additional N-terminal fragment. Furthermore, a non-structural protein (NS-1) is synthesized, which is required for virus replication.

Clinical spectrum of a Parvovirus B19 infection:

- Fifth disease (erythema infectiosum): It progresses in 20% of cases without symptoms or with flu-like symptoms. Apart from the exanthema, polyarthralgia and a generalized swelling of the lymph nodes can occur occasionally.
- Aplastic crisis in patients suffering from hemolytic anemias.
- Spontaneous abortion/hydrops fetalis after infection during pregnancy.
- Persistent infection in immunosuppressed persons with chronic anemia and bone marrow depression.

By the kind of antigen pattern and the avidity determination of IgG antibodies the *recomLine* Parvovirus B19 provides an indication of the time of infection. The determination of antibodies against the NS-1 antigen can be helpful in the clarification of persistent Parvovirus B19 infections. With the use of VP-2 particles the presentation of conformational epitopes in addition to the linear epitopes is achieved.




Product advantages

- Recombinant antigens
 - High sensitivity and specificity
 - Easy and clear interpretation due to easy to read bands
- Easy test procedure; automation possible
- Easy and objective strip interpretation and documentation with *recomScan* software
- Test procedure and reagents identical in all MIKROGEN strip tests - reagents exchangeable
- Safe evaluation due to strip specific controls (cut-off and conjugate control)
- Separate detection of IgG and IgM antibodies possible
- Easy and reliable determination of avidity possible
- CE label: The *recomLine* Parvovirus B19 tests meet the high standard of the EC directive 98/79/EC on in vitro diagnostic medical devices
- **Antigen reaction pattern provides accurate determination of infection status**

Related recombinant Parvovirus B19 antigens

| Parvovirus antigen | Recombinant antigen | Size [kDa] |
|--|---------------------|------------|
| Main capsid antigen (conformational epitopes) | VP-2p | 56 |
| N-terminal parts of structure proteins VP-1 and VP-2 | VP-N | 60 |
| VP-1 specific segment (differentiation to VP-2) | VP-1S | 31 |
| Main capsid antigen (linear epitopes) | VP-2r | 56 |
| C-terminal part of structure proteins VP-1 and VP-2 | VP-C | 42 |
| None structural protein NS-1 | NS-1 | 75 |

Test Principle and Procedure



- 1st Incubation** A test strip loaded with Parvovirus B19 antigens is incubated with diluted serum or plasma in a dish for **1 hour**.
- wash 3 times
- 2nd Incubation** Peroxidase conjugated anti-human antibodies (IgG or IgM specific) are added. Incubate for **45 minutes**.
- wash 3 times
- Color reaction** **8 minutes** after addition of the coloring solution, insoluble colored bands develop at the sites on the test strips occupied by antibodies.

Examples

| date | day after symptoms | recomWell Parvovirus B19 | | recomLine Parvovirus B19 | | | | | | | strip result | interpretation | |
|---------------|--------------------|--------------------------|-------------------|--------------------------|-------|------|-------|-------|------|------|--------------|----------------|------------------------|
| | | IgG | IgM | Cutoff | VP-2p | VP-N | VP-1S | VP-2r | VP-C | NS-1 | | | |
| 1 | 2 | 164 U/ml positive | 102 U/ml positive | IgG | ++ | + | - | + | - | | | positive | acute infection |
| | | | | avidity | n | n | | | | | | low avidity | |
| | | | | IgM | ++ | ++ | +/- | ++ | - | | | positive | |
| 2 | 17 | 251 U/ml positive | 84 U/ml positive | IgG | ++ | +++ | ++ | +++ | ++ | | | positive | acute infection |
| | | | | avidity | i | n | | | | | | low avidity | |
| | | | | IgM | ++ | ++ | +/- | ++ | - | | | positive | |
| 3 | 59 | 296 U/ml positive | 33 U/ml positive | IgG | ++ | +++ | +++ | ++ | + | | | positive | status after infection |
| | | | | avidity | h | h | | | | | | high avidity | |
| | | | | IgM | + | + | - | - | - | | | positive | |
| 4 | 144 | 219 U/ml positive | 21 U/ml indeterm. | IgG | ++ | +++ | +++ | ++ | - | | | positive | past infection |
| | | | | avidity | h | h | | | | | | high avidity | |
| | | | | IgM | +/- | - | - | - | - | | | negative | |
| blood donor 1 | | positive | negative | IgG | +++ | +++ | +++ | +/- | - | | | positive | past infection |
| | | | | avidity | h | h | | | | | | high avidity | |
| blood donor 2 | | positive | negative | IgG | ++ | - | - | - | - | | | positive | past infection |
| | | | | avidity | | | | | | | | no conclusion | |

Typical Parvovirus B19 course of infection and two examples of healthy blood donors with past infections

The course of infection covers a period of approx. 21 weeks after appearance of symptoms. IgM reactivity decreased significantly between sample 2 and 3. The decrease of the VP-C IgG reactivity, the high avidity and the appearance of NS-1 IgG antibodies indicate a status after infection (sample 3). High avidity appears earliest 4 weeks after infection. NS-1 antibodies titers appear only six to eight weeks after infection, but not in all patients. In a study samples from blood donors were tested, 66% were IgG positive. Many samples presented this typical antigens pattern in past infections as shown. In Germany the infection rate in adults is about 60% and in children between 15-30%.

Article-No

4472 **recomLine Parvovirus B19 IgG [Avidity]***
Reagents for 20 determinations

4473 **recomLine Parvovirus B19 IgM**
Reagents for 20 determinations

11010 **Avidity reagent**
Reagent for 25 avidity determinations

*[] optional available as additional reagent

Storage

At +2°C - +8°C