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
## Instructions for use

### Composition

Pathogen-specific positive material in protein containing phosphate buffer with Tween 20; liquid, 1 x 3 ml; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus-surface antigen) and anti-HCV Ab; preservative: < 0.1 % sodium azide; colouring: 0.01 g/l Bromphenol blue.

### Intended Use

SERION ELISA control are positive controls for monitoring qualitative and quantitative antibody determination when using immunoassays from the SERION ELISA product line. They are to be regarded as additional controls to the reagents supplied with the test kits. They are used to determine validity of SERION ELISA test runs as well as precision and reliability of the method. SERION ELISA control are particularly recommended as an aid for internal quality management in accredited laboratories.

### Stability and Storage

The ready-to-use SERION ELISA control are stable for 24 months from date of production when stored at 2–8 °C; expiry date see label. Do not freeze.

### Precautions

Although all SERION ELISA control have been found to be negative for HBsAg, HCV and HIV-antibodies, they should be considered potentially infectious. The product contains sodium azide as a preservative. Swallowing and contact with skin or mucous membranes must be avoided!

### Target Values and Target Ranges

In comprehensive validation studies, target values and target ranges for each test-specific SERION ELISA control are determined and documented on lot-specific certificates of analysis. These are recommended as reference values for the quality management in accredited laboratories when using SERION ELISA control with the corresponding SERION ELISA kits. SERION ELISA control generally show positive results in test runs with SERION ELISA. However, a few SERION ELISA tests were evaluated to correspond to so-called clinical threshold, e.g. in order to exclude high seroprevalences. Therefore, the evaluation of individual SERION ELISA control might deliver borderline results in exceptional cases.

### Preparation

SERION ELISA control are ready-to-use and must not be diluted any further. Rf-absorption must not be performed. SERION ELISA control must be brought to room temperature before using. Before pipetting, SERION ELISA control should be agitated (e.g. by using a vortex mixer) in order to ensure a homogenous solution.

### Performance

In course of the test run pipette 100 µl of SERION ELISA control into the cavity of the corresponding SERION ELISA test strip. SERION ELISA control are recommended to be used in each SERION ELISA test run and the results recorded on a control chart. SERION ELISA control are recommended for use with the corresponding SERION ELISA test. Please also refer to the instruction manual of the SERION ELISA.

### Limitations

SERION ELISA control must not be used after expiry date. In case of microbiological contamination or turbidity, the product should be discarded. SERION ELISA control are independent from the controls included in the SERION ELISA kit. They are not suitable for use as standard. SERION ELISA control can be used with other tests only after validation by the user.

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